



# GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. – Via Marconi, 1 - 20060 Gessate (MI) Italy  
Italia: tel 199 400 401 – fax 199 400 403 – [gima@gimaitaly.com](mailto:gima@gimaitaly.com)  
Export: phone + 39 02 953854209/221/225 – fax + 39 02 95380056  
- [export@gimaitaly.com](mailto:export@gimaitaly.com) – [www.gimaitaly.com](http://www.gimaitaly.com)

## Code 27258

### Electronic medium pet scale

**Thank you for choosing our product.**  
**Your new electronic petscale was designed for measuring your pet's weight for a many-year-long time period.**

**Please, read the instruction for use before using the scale.**

#### Technical specification:

1" LCD display  
Capacity: 100kg /220Lbs  
Graduation: 100g / 0,2Lb  
Auto off  
Auto overload indicator ( "ERR" )  
Low battery indicator ( "LO" )  
kg, lb, St unit of measurement can be switched

#### User Guide

1. Before using the scale for the first time, please remove isolating tab from the battery at the backside of the scale.
2. Place the scale on a flat, hard level surface. On a soft surface ( for example: on rug ) the scale may show inaccurate value.
3. To start the scale, tap the scale platform with your foot.
4. For the first time the display shows a running "8" and after "0,0" appears on the display.
5. Wait until "0,0" is displayed!
6. Place pet on the scale.
7. Then your pet's weight will appear on the display flashing two times. The measured weight can be seen for about 10 seconds.
8. After weighing, it will automatically turn off approximately in 6 seconds.

#### Power:

The scale run on replaceable Lithium battery. ( CR 2032 ) included

#### Replacing batteries:

When "LO" appears on the display, replace the batteries.

#### Error display:

Weighing above 100 kg, "ERR" is displayed.  
Please do not overload the scale! In case of overloading the scale get damaged. In this case we have no liability for the product.

#### Environment:

The environmental temperature may influence the weighing accuracy. You can get correct weighing accuracy at room temperature (20-25°C)

**Remarks:**

Please do not use the scale in the immediate vicinity of strong emitters such as mobile phones or microwave ovens. To clean the scale use damp clothes without any solvent.



Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

**GIMA WARRANTY CONDITIONS**

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA. During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the warranty. The repair or replacement performed during the warranty period shall not extend the warranty. The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use. GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc. The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed. The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.