

MONOPOLAR ELECTROSURGERY UNIT - 50W



• 30516 DIATERMO 106 - 50W - monopolar

This HF electrosurgery unit, is particularly suitable for small operations and is recommended for dentistry, first aid and dermatology work. It offers precision of cut and immediate coagulation and give excellent results in a large range of applications: from uterine coagulation to microcoagulation in dermatology. The modes of operation, CUT, COAG and



TECHNICAL SPECIFICATIONS

Output power:	CUT	50W/400 ohm	
	COAG	40W/400 ohm	
	BLEND	45W/400 ohm	
	MICRO	-	
Working frequency:	600 KHz +/- 10%		
Size:	180x180x h 100 mm		Weight: 2.8 kg
Operating Voltage:	115/230 V - 50/60 Hz switchable		
Electric input power:	90 VA		Fuse: 2x1 AT
Alarm for overheating		Footswitch operating	
Norms:	EN60601-1 Class I CF, EN 60601-1-2 93/42/EEC (Class IIb)		
Made in Italy		Warranty: 12 months	

STANDARD ACCESSORIES

KIRKY single-use handle without buttons
 Standard foot-switch (30568)
 Rod neutral electrode and cable
 Kit of 6 short autoclavable electrodes, Ø 2.4 mm, length 5 cm
 Power cable 2 m - Shuko plug
 User manual: available in English, Italian, Spanish. French



BLEND, can be easily selected by a switch on the front panel. The complete isolation between patient and the unit and strict respect of international norms guarantee total safety.