



# SU-180 Urine Analyzer

Urine Analyzer

# PROKAN



# SU-180 Urine Analyzer

## SPECIFICATIONS

- Large LCD Display, Navigation System, Chinese and English Operation interface, and data outputting
- Cold light source with high luminosity, good stability, long shelf life
- Innovative light sensor, high sensitivity, strong anti-Jamming
- Automatic calibrate the impact of uric acid basicity rate to the test result
- Fast speed 120 test per hour
- Single, Continuous mode selection
- Standard RS-232 interface for connection with LIS or urine sediments analyzer
- Storage of 1000 test results
- Automatic exception setting
- Built-in thermal printer
- Green designing test strip cutting block for easy cleaning



## MAIN TECHNICAL SPECIFICATIONS

- Test item: GLU, BIL, SG, PH, KET, BLD, PRO, URO, NIT, LEU, Vc
- Wave length of monochromatic light: 470nm, 540nm, 620nm
- Principle: Reflectance photometer
- Strips: 11 items.
- Throughput: 120 tests/hour, 60 test /hour and manual trigger
- Display: Large LCD display, Resolution: 240\*128
- Memory: 1000 results
- Printer: Built in thermal printer, paper width: 57mm
- Interface: Standard RS-232
- Working environment: Temperature 15°C - 30°C Humidity < 70%
- Power Supply: 100V-240V 50~60HZ

PROKAN ELECTRONICS INC.

 **ALL4VETS**  
PRE VETERINÁRNÝCH LEKÁROV

Výhradný distribútor pre Česko a Slovenskú republiku:  
All4vets spol. s r.o., Potočná 14, 921 01 Piešťany  
www.all4vets.eu, e-mail: info@all4vets.eu  
Tel.: +421 915 835 345, +421 915 974 046