


SN-T2/SN-T2R

Syringe Pump



Convenient: portable

Intelligent: Syringe calibration
Wireless communication

Appearance :500g only

Safety: Two sensors: Pressure
sensor, Position sensor



All4vets, s.r.o.
Radlinského 37, 921 01 Piešťany
info@all4vets.eu

SN-T2/SN-T2R

Specification:

Model	
Flow rate	0.1 ml/h ~ 1800 ml/h
Accuracy of infusion	±2% ~ ±5%
Bolus rate	5.0 ml/h ~ 1800 ml/h
Delivered volume	0.1 ml ~ 9999 ml
History records	1000PCS
Delivery limit	0.1 ml ~ 9999 ml
KVO rate	0.1 ml/h ~ 5.0 ml/h
Syringe	5ml , 10ml , 20ml , 30ml , 50/60ml
Syringe Specification	30 brands of syringe
Occlusion detection pressure	H,C,L
Alarm sound level	Adjustable---3 levels
Button sound	Adjustable---3 levels, can be on/off
Working mode	Speed mode, Time mode, Body weight mode, Infusion program
Other function	PCA(SN-T2R), Display of remained volume
Communication	Wireless with computer (SN-T2R)
Power	Support more than 10 hours at the flow rate 5ml/h
Power supply	DC 12V
Battery	2X9V dry battery
Voltage	12VA
Working condition	Relative humidity: 20%~80%; Ambient temperature: +5°C~+40°C ; Ambient pressure: 76kPa ~ 106kPa.
Weight (with clamp)	0.5 kg

Function:

Dual-CPU monitoring function	Automatic slow impact dose release after occlusion	KVO function	Abundant optic and sound alarm
Preset volume	Display of remained volume	Fast rate control	External DC power interfere
Syngce calibration	Key-lock function	Preset volume	Wireless communication function(SN-T2R)
Automatic identification of syringe specifications	Delivered volume	Independently PCA button with extension wire	
3-level occlusion alarm limits setting	1000 history records		

Predaj a distribúcia:

All4vets, s.r.o.

Radlinského 37, 921 01 Piešťany

info@all4vets.eu

+421 907 145 915

