

Portable Slit Lamp Microscope

SLM-6L/6M/6H Portable Slit Lamp Microscope

➤ Configuration(SLM-6L) Portable Slit Lamp

with 10X eyepieces

with Li-ion Battery, Battery Charger

Cena přístroja : 1470 € bez DPH

With metal box

With focus stick, Forehead Bracket

➤ Configuration(SLM-6M) Portable Slit Lamp

Changing magnification by rotating objective lens (patent)

with 10X eyepieces

with Li-ion Battery, Battery Charger

Cena přístroja : 1585 € bez DPH

With metal box

With focus stick, Forehead Bracket

➤ Configuration(SLM-6H) Portable Slit Lamp

with 10X,16X eyepieces

with Li-ion Battery, Battery Charger

Cena přístroja : 1693 € bez DPH

With metal box

With focus stick, Forehead Bracket



Parameters:

Model		SLM-6L	SLM-6M	SLM-6H
Microscope	Type	Converging		
	Objective Lens	10X	10X, 16X	10X, 16X
	Eyepiece	10X	10X	10X, 16X
	Total Magnifications	10X	10X, 16X	10X, 16X, 25X
	Working Distance	81mm		
	Pupil Distance	45mm~78mm		
	Diopter Compensation	-8D~+8D		
Illumination	Slit Width	0-12mm Continuous		
	Slit Length	12mm	1-12mm	1-12mm
	Slit Angle	Horizontal +30°		
	Filter	Heat absorption, yellow, red free(green), Cobalt blue		
	Maximum Illumination	>30000 Lux		
	Power	5VA		
Power Source	LED working Voltage	DC3V, 350mAh		
	Battery	Li-ion battery DC7.4V/680mAh, Full battery can supply continuous working time 5-6 hours		

Portable Slit Lamp Microscope

YZ-3 Portable Slit Lamp Microscope

Features:

- full-featured, lightweight and durable and is suitable for mobile medical, community hospitals and wards checkup
- Excellent optical performance, total 800g weight, easy to operate and practical enough
- Converging microscope provides sharp observation.

Cena přístroja : 1756 € bez DPH



Parameters:

Microscope Type	Converging microscope
Objective Magnifications	Two steps magnifications
Magnifications	10X,16X
Diopter Adjustment	±8D
Pupil Distance	50mm~75mm
Working Distance	60mm
Slit Width	0~12mm continuously adjustable
Diameter of Light Spot	0.2mm,1mm,5mm,12mm
Filter	Color Temperature Compensation, Red free, Cobalt Blue
Slit Angle	±30°
Light Source	LED
Maximal Illumination	15000Lx
Power	7.4V 2200mA Li Battery, Rechargeable
Working Time	2.5 hours for common use
Packing Volume	216mm X 100mm X 225mm
Net Weight	800g

Portable Slit Lamp Microscope

SLM-1A/B Portable Slit Lamp Microscope

Features:

- Advanced optics, 10X& 16X magnification
- Controllable illumination from maximum to zero
- Enhanced portability during examination

Cena přístroja SLM-1A : 2450 € bez DPH

Cena přístroja SLM-1B: 2540€ bez DPH



Parameters:

Model		SLM-1(A)	SLM-1(B)
Microscope	Type	Binocular Portable Slit Lamp	
	Magnification	10X	10X, 16X
	Eyepiece	Convergence binoculars, 13°	
	Working Distance	80mm	
	Pupil Distance	45mm~78mm	
	Diopter Compensation	-8D~+8D	
Illumination	Slit Width	0.15mm, 0.5mm, 0.8mm, 1.6mm, 12mm, 1X1mm Block	
	Slit Length	12mm	
	Slit Angle	Horizontal +-60°	
	Filter	Heat absorption, Color Temperature Compensation, red free, Cobalt blue	
	Maximum Illumination	>=50000 Lux	
	Light Source	White LED	
Power Source	LED working Voltage	DC7.4V, 950mAh	
	Battery	Li-ion battery DC7.4V/950mAh, continuous working time>=3hours after fully charged	

Portable Slit Lamp Microscope

KJ5S1 Portable Slit Lamp Microscope

Features:

- Portable slit lamp, one hand operated
- With all Functions of desktop slit lamp but ideal for pediatrics, Seniors, nursing home, carry-out or mobile hospital.
- This handheld slit lamp has four light spots: 1mm, 3mm, 5mm, 10mm
- Fully charged can be used for about 6 hours
- Smart design, Net Weight 750g, easy carry with

Cena přístroja : 1712 € bez DPH



Parameters:

Microscope type:	converging microscope
Objective magnification:	1x
Eyepiece magnification:	10x, 16x
Diopter compensation	+/-7D
Pupil distance:	80mm
Slit width:	0-10mm continuously adjustable
Diameter of light spot:	1mm, 3mm, 5mm, 10mm
Filter:	Heat-absorption, Grey, Red-free, Cobalt blue
Slit angle:	horizontal +30
Light source:	white LED
Power:	7.2V 750Ah Li Battery, Rechargeable
Working Time:	5-6 hours
Packing Volume:	40cm*20cm*23cm
Net weight:	750g
Gross weight:	2600g

Portable Slit Lamp Microscope

KJ5S2 Portable Slit Lamp Microscope

Features:

- Portable, one hand operated
- With all Functions of desktop slit lamp but ideal for pediatrics, Seniors, nursing home, carry-out or mobile hospital.
- light spots: 10mm
- Fully charged can be used for about 6 hours
- Smart design, Net Weight 750g, easy to carry with
- 16 months warranty
- CE, ISO, FDA



Cena prístroja : 1510 € bez DPH

Parameters:

Microscope Part:	Microscope Type:	Biocular Crossed Stereomicroscope
	Objective Magnification:	1X
	Eye piece Magnification:	10x, 16x
	Dia of Light Spot (mm):	10
	Pupil Distance (mm):	48-70
	Adjustable Diopter:	-7D+7D
	Working distance:	80mm
illumination:	Slit width:	0~10mm adjustable
	Slit height:	1~10mm adjustable
	Slit Angle :	horizontal +/- 30°
	Filter:	Heat-absorption, Grey, Red free, Cobalt blue
	Light source:	High illuminance white LED
	Working Time:	5-6 hours
Power:	7.2V Li-ion battery 750 Ah, Rechargeable	
Charger:	Input voltage:	AC 100-240V 60-50Hz
	Output voltage:	DC 18V 500mAh

Portable Slit Lamp Microscope

KJ5S3 Portable Slit Lamp Microscope

Features:

- New developed digital portable slit lamp is very convenient to take photo from patients.
- Can be equipped with Nikon camera
- Portable, one hand operated
- With all Functions of desktop slit lamp but ideal for pediatrics, Seniors, nursing home, carry-out or mobile hospital.
- Multi light spots: 1mm, 3mm, 5mm, 10mm
- 16 months warranty
- CE, ISO



Cena přístroja : 2265 € bez DPH

Parameters:

Microscope Part:	Microscope Type:	Biocular Crossed Stereo-microscope
	Objective Magnification:	1X
	Eye piece Magnification:	10x, 16x
	Dia of Light Spot(mm):	1, 3, 5, 10
	Pupil Distance (mm):	48-70
	Adjustable Diopter:	-7D+7D
	Working distance:	80mm
illumination:	Slit width:	0~10mm adjustable
	Slit height:	1~10mm adjustable
	Slit Angle :	horizontal+/- 30°
	Filter:	Heat-absorption, Grey, Red free, Cobalt blue
	Light source:	High illuminance white LED
	Working Time:	5-6 hours
Power:	7.2V Li-ion battery 750 Ah, Rechargeable	
Charger:	Input voltage:	AC 100-240V 60-50Hz
	Output voltage:	DC 18V 500mAh

Portable Slit Lamp Microscope

ML-1 Portable Slit Lamp Microscope

Features:

- Portable Pocket Optical Hand-held Slit Lamp Biomicroscope for Anterior Eye Examination
- Lightweight, single hand operate;
- Link-able Eyepieces;
- Cordless and rechargeable, no Restriction of Time and Space for Eye Examination
- Packed in Alu-Alloy Carrying Case, with black and white colors;
- It can be digitized by Coopix Camera or Phone Eyepiece Adapters of iPhone, Samsung, iPad, iPod etc.



Cena prístroja : 1330 € bez DPH

Parameters:

Microscope Part	
Eyepiece magnification:	10x,16x
Objective Lens	1x
Total Magnifications	10X,16X
Diopter compensation	+/-5D
Pupil distance:	50~72mm
Working Distance	80mm
Slit Illumination Part	
Slit width:	0-10mm continuously adjustable
Slit height	0.5mm,3mm,5mm,10mm
Diameter of light spot:	0.5mm,3mm,5mm,10mm
Filter:	Heat absorption, Grey, Red-free, Cobalt blue
Slit angle:	horizontal +/-30°
Light source:	white LED
Power:	7.2V 750Ah Li Battery,Rechargeable
Working Time:	5-6 hours
Voltage/Frequency	AC 110V+/-10% , 60Hz 1Hz AC 220V+/-10% , 50Hz 1Hz
General	
Packing Volume:	41cm*20cm*23cm
Net weight:	750g
Gross weight:	2.6kg